




Multitrim User Manual

Kenji Gerhardt

Multitrim Workflow



Requirements and Installation



Multitrim requires the python programming language and the following tools to be installed:

- FaQCs:
 - <https://github.com/LANL-Bioinformatics/FaQCs>
- Fastp:
 - <https://github.com/OpenGene/fastp>
- Falco:
 - <https://github.com/smithlabcode/falco>
- SeqTK:
 - <https://github.com/lh3/seqtk>

The easiest way to do all of this is to use Conda. Instructions for doing so are on the next few slides.

Installation of Miniconda



1. (Windows only): Get [Ubuntu for Windows](#) and install it.
2. Download the appropriate version of [Miniconda](#):
3. Follow the installation instructions:
 - a. [Ubuntu for Windows/Linux](#)
 - b. [MacOS](#)

Installation of Multitrim



1. Acquire multitrim with git by running the following command in your terminal:
 - a. `git clone https://github.com/KGerhardt/multitrim`
2. Create the multitrim environment by running the following command:
 - a. `conda env create -f multitrim/multitrim.yml`
3. Activate the conda environment for multitrim with the following command:
 - a. `conda activate multitrim`
4. You're ready to run the multitrim script.
 - a. You need to rerun (3) each time you open your terminal. Conda environments are not active by default.

Quality of Life (Optional)



- Make multitrим available as an alias by running the following command in your terminal:
 - `alias multitrим="python path/to/multitrим_directory/multitrим/multitrим.py"`
- Make it permanent by adding it to your profile:
 - [Ubuntu for Windows/Linux:](#)
 - [MacOS](#)

If you don't do this, you'll need to remember where the multitrим.py script is located and use it as `python /path/to/multitrим_directory/multitrим/multitrим.py [OPTS]` each time. It's much easier to create the alias once and then always have it available when you open your terminal.

Using Multitrim - Basic usage



Note: I will assume going forward that you've followed the optional instructions on slide 5, and can use multitrim by simply typing “multitrim” on the command line.

- Single end reads:
 - `multitrim -u [unpaired_reads.fq]`
- Paired end reads:
 - `multitrim -1 [forward_reads.fq] -2 [reverse_reads.fq]`

That's it.

Note: input Reads can be gzipped (they probably end with “.fq.gz”), and multitrim still works.

Using Multitrim - Controlling Outputs



- Multitrim has a standard set of naming conventions for its outputs:
 - (directory)/(prefix_)(file_identifier)_(base_name).(extension)
- There is no default prefix, and the default directory is the current working directory.
 - You can specify a directory with the -d option
 - You can specify a prefix with the -p option
- The naming convention for reads and their QC reports include information from the input file names.
 - This means you can run multitrim on multiple sets of reads, place the results in the same directory, and all critical information will be kept.
 - If you want to keep the fastp reports, the detected adapter files, and the subsample results for multiple reads in a single directory, use prefixes to make these names unique.

Using Multitrim - Trimming Options



Multitrim uses both FaQCs and fastp to trim. FaQCs performs an initial trim, and then fastp trims those results. The following settings for these tools are available to a user of multitrim:

- FaQCs:
 - `--score [INT]`; sets FaQCs score parameter. Default 27.
- Fastp:
 - `--advanced`; applies poly-G tail trim and low-complexity read filters. Read fastp documentation for more detail
 - `--window [INT]`; sets the width of the sliding window for fastp trimming. Default 3 bp.
 - `--window_qual [INT]`; sets the minimum avg. quality for bases in the window. Default 20.
 - `--min_L [INT]`; minimum read length post-trim for a read to be kept. Default 50.

Using Multitrim - Other Useful Options



- `--threads [INT]`
 - Uses the specified number of threads to perform trimming and QC. More threads = faster. Default 1.
- `--max`
 - Uses all of the cores available to your computer/HPC job. Overwrites `--threads` if both are supplied.
 - **You should always use this option.**
- `--min_adapt_pres [float]`
 - Adapters must be detected in the specified percent of reads during subsampling to be considered present. Default 0.1

Using Multitrim - Other Less-Useful Options



- `--phred_fmt`
 - Exists to resolve an issue with FaQCs in SE mode. Sets phred scoring offset. Default 33.
- `--skip_faqcs`
 - Skips the FaQCs trimming step. Generally results in over-trimmed reads.
- `--skip_fastp`
 - Skips the Fastp trimming step. Generally results in under-trimmed reads.
- The following allow you to supply the location of binaries for each tool, in case you cannot conda install
 - `--falco`
 - `--seqtk`
 - `--faqcs`
 - `--fastp`

Multitrim Outputs:



- Both SE and PE:
 - detected_adapters.fasta
 - post_trim_fastp.html
 - post_trim_fastp.json
 - Subsample_Adapter_Detection.stats.txt
 - Subsample_Adapter_Detection_qc_report.pdf
- SE only:
 - unpaired.pre_trim_qc_[file_base_name].html
 - unpaired.post_trim_qc_[file_base_name].html
 - unpaired.trimmed_[file_base_name].fq.gz
- PE only:
 - 1.pre_trim_qc_[forward_base_name].html
 - 1.post_trim_qc_[forward_base_name].html
 - 1.trimmed_[forward_base_name].fq.gz
 - 2.pre_trim_qc_[reverse_base_name].html
 - 2.post_trim_qc_[reverse_base_name].html
 - 2.trimmed_[reverse_base_name].fq.gz

Note: all files here will begin with the user's chosen prefix, if one is supplied with -p

Note: File base names are created as such:
base_name.extension(.gz) -> base_name

Understanding the outputs:



- `*pre_trim_QC*`
 - Falco QC report on your inputs before any trimming
- `*post_trim_QC*`
 - Falco QC report on the final, fully trimmed outputs
- `[1/2/unpaired].trimmed*.fq.gz`
 - Your trimmed reads as a gzipped FASTQ file
- `detected_adapters.fasta` - FASTA file containing detected adapter sequences.
 - May be empty if no adapters are found.
- `post_trim_fastp.[html/json]`
 - graphical/text report from fastp trim.
- `Subsample_Adapter_Detection.[pdf/stats.txt]`
 - graphical/text reports from FaQCs on the subsampled 100K reads used to detect adapters.